


## Face Mask List


(1 / 4)

### ■ Plane Type

Image	No	Model	Useage	Description
	1	DM-ES-A-W	Family Mask	Outer layer cloth color : White Used for Adults
	2	DM-ES-A-B	Family Mask	Outer layer cloth color : Blue Used for Adults
	3	DM-ES-A-G	Family Mask	Outer layer cloth color : Green Used for Adults
	4	DM-ES-A-P	Family Mask	Outer layer cloth color : Pink Used for Adults
	5	DM-ES-A-Y	Family Mask	Outer layer cloth color : Yellow Used for Adults
	6	DM-ES-K-W	Family Mask	Outer layer cloth color : White Used for Kids
	7	DM-ES-K-B	Family Mask	Outer layer cloth color : Blue Used for Kids
	8	DM-ES-K-G	Family Mask	Outer layer cloth color : Green Used for Kids
	9	DM-ES-K-P	Family Mask	Outer layer cloth color : Pink Used for Kids
	10	DM-ES-K-Y	Family Mask	Outer layer cloth color : Yellow Used for Kids

※ These are ear suspension type.

### ■ C type


Image	No	Model	Useage	Description
	1	AP069-PS1	Medical Mask EN149 FFP1	Particle Filtration Efficiency (PFE) > 90% Bacterial Filtration Efficiency (BFE) > 99.9% Class : CE EN149:2001 FFP1NR Hands over the parallel suspension type
	2	AP089-PS2	Medical Mask EN149 FFP2	Particle Filtration Efficiency (PFE) > 95% Bacterial Filtration Efficiency (BFE) > 99.9% Class : CE EN149:2001 FFP2NR Hands over the parallel suspension type
	3	AP708-PS3	Medical Mask NIOSH N95	Particle Filtration Efficiency (PFE) > 95% Bacterial Filtration Efficiency (BFE) > 99.9% Class : NIOSH N95 TC-84A-4990 Hands over the parallel suspension type

※ The color of these masks is white.


## Face Mask List

( 2 / 4 )

### ■ 3D-JP Type

Image	No	Model	Useage	Description
	1	3DMJP-ES-A	Family Mask Medical Mask	Particle Filtration Efficiency (PFE) > 95% Bacterial Filtration Efficiency (BFE) > 99.9% Class : NIOSH N95 Outer layer cloth color : 4 colors ( White · Blue · Green · Pink ) Used for Adults Ear suspension type
	2	3DMJP-ES-K	Family Mask Medical Mask	Particle Filtration Efficiency (PFE) > 95% Bacterial Filtration Efficiency (BFE) > 99.9% Class : NIOSH N95 Outer layer cloth color : 4 colors ( White · Blue · Green · Pink ) Used for Kids Ear suspension type

### ■ 3D-GL Type

Image	No	Model	Useage	Description
	1	3DMGL-ES-A-W	Family Mask	Outer layer cloth color : White Used for Adults
	2	3DMGL-ES-A-B	Family Mask	Outer layer cloth color : Blue Used for Adults
	3	3DMGL-ES-A-P	Family Mask	Outer layer cloth color : Pink Used for Adults
	4	3DMGL-ES-A-Y	Family Mask	Outer layer cloth color : Yellow Used for Adults
	5	3DMGL-ES-K-W	Family Mask	Outer layer cloth color : White Used for Kids
	6	3DMGL-ES-K-B	Family Mask	Outer layer cloth color : Blue Used for Kids
	7	3DMGL-ES-K-P	Family Mask	Outer layer cloth color : Pink Used for Kids
	8	3DMGL-ES-K-Y	Family Mask	Outer layer cloth color : Yellow Used for Kids

※ These are Ear suspension type

## Face Mask List

( 3 / 4 )

### ■ Cone type

Image	No	Model	Useage	Description
	1	E520-PS-N95	Work mask NIOSH N95	Particle Filtration Efficiency (PFE) > 95% Outer layer cloth color : White Class : NIOSH N95 TC-84A-4362 Used for Adults Hands over the parallel suspension type
	2	E220V-PS-W2	Work mask EN149 FFP2	Particle Filtration Efficiency (PFE) > 95% Outer layer cloth color : White With valve Class : EN149: 2001-FFP2 Used for Adults Hands over the parallel suspension type
	3	F550-PS-N95	Work mask NIOSH N95	Particle Filtration Efficiency (PFE) > 95% Outer layer cloth color : White Class : NIOSH N95 TC-84A-4394 Used for Adults Hands over the parallel suspension type
	4	AP810-PS-N95	Work mask NIOSH N95	Particle Filtration Efficiency (PFE) > 95% Outer layer cloth color : White Class : NIOSH N95 TC-84A-5496 Used for Adults Hands over the parallel suspension type
	5	AP088-PS-W2	Work mask EN149 FFP2	Particle Filtration Efficiency (PFE) > 95% Outer layer cloth color : White Class : CE EN149:2001 FFP2NR Used for Adults Hands over the parallel suspension type

## Face Mask List

( 4 / 4 )

### ■ Nano Antibacterial Face Mask

Image	No	Model	Useage	Description
	1	<b>NAFM-CFS-A</b>	Family Mask Medical Mask Work mask	Particulate Filtration Efficiency (PFE)Ave.99.5% Bacterial Filtration Efficiency (BFE)> 99.9% <b>Approved : PM2.5 A Grade</b> Protection Efficiency of Mask : $\geq 90\%$ Sub-Micron Filtration Efficiency (NaCl) : $\geq 95\%$ PM2.5 Filtration Efficiency (KCI) : $\geq 99\%$ PM2.5 Filtration Efficiency (PAO) : $\geq 99\%$ With PP Non-Woven Sputtered by Ag/Cu/Ti Color : Outer layer cloth color : Green Ear Suspension with Back Hook
	2	<b>NAFM-CFS-B</b>	Family Mask Medical Mask Work mask	<b>Approved : PM2.5 B Grade</b> Protection Efficiency of Mask : $\geq 85\%$ Sub-Micron Filtration Efficiency (NaCl) : $\geq 90\%$ PM2.5 Filtration Efficiency (KCI) : $\geq 95\%$ PM2.5 Filtration Efficiency (PAO) : $\geq 95\%$ With PP Non-Woven Sputtered by Ag/Cu/Ti Color : Outer layer cloth color : Green Ear Suspension with Back Hook
	3	<b>NAFM-CFS-C</b>	Family Mask Medical Mask Work mask	<b>Approved : PM2.5 C Grade</b> Protection Efficiency of Mask : $\geq 75\%$ Sub-Micron Filtration Efficiency (NaCl) : $\geq 80\%$ PM2.5 Filtration Efficiency (KCI) : $\geq 90\%$ PM2.5 Filtration Efficiency (PAO) : $\geq 75\%$ With PP Non-Woven Sputtered by Ag/Cu/Ti Color : Outer layer cloth color : Green Ear Suspension with Back Hook
	4	<b>NAFM-CFS-D</b>	Family Mask Medical Mask Work mask	<b>Approved : PM2.5 D Grade</b> Protection Efficiency of Mask : $\geq 65\%$ Sub-Micron Filtration Efficiency (NaCl) : $\geq 80\%$ PM2.5 Filtration Efficiency (KCI) : $\geq 90\%$ PM2.5 Filtration Efficiency (PAO) : $\geq 75\%$ With PP Non-Woven Sputtered by Ag/Cu/Ti Color : Outer layer cloth color : Green Ear Suspension with Back Hook

\* These are used for Adults.